


# BioQube



Revolutionizing blood testing with  
Quantum Diamonds



# Our Vision

- **BioQube** seeks to provide a new and exciting future for biodiagnostics. Our solution combines ease of use with pinpoint precision without requiring any invasive procedures.
- **BioQube's** Vision is to solve problems in the diagnostics of almost any disease. Alzheimer is the first step on this mission.

# Problem Space

## Speed & Usability

Current Tests are often slow and require highly trained staff, which is not readily available everywhere.


## Cost & Flexibility

Due to high manual labor usage, costs of current tests are high and they can not be carried out everywhere on a flexible timescale.

## Accuracy

Current tests are often not accurate enough to diagnose certain diseases early enough for effective treatment.

# Our **Solution** powered by Quantum NV Diagnostics



FASTER

With our NV-Diamond setup utilizing microfluidic Chips, we can provide fast testing, with results only limited by the computational power needed for processing of the measurements.



COST EFFECTIVE

BioQube requires no skilled labor. Additionally, due to the reusable nature of our Design, we have no single use components. This leads to repeat testing being very cost effective.



MORE  
SENSITIVE


BioQube leverages Quantum NV-Vacancy technology to achieve unprecedented sensitivity and accuracy.

## The **BioQube**

Enabling easy fast, low cost and accurate blood diagnosis without the need for specialised staff.




# Our **Solution** the Sample Handling



The Microfluidic Chip provides easy, fully automated and reliable sample preparation and handling.

The BioQube folds open giving easy access for up to 3 different blood samples to be placed into the machine. Then the sample and the nanoparticles needed for the testing are fed through and mixed in the Tesla mixer. The mixed and prepared sample is then passed on to the NV-Chip for the diagnostics. The Chip has 7 compartments for different nanoparticles so each sample can be tested for up to 7 different diseases per diagnostic run. After the analysis has been done the Chip flushes itself with a cleaning solution, making it ready for the next run.


# Our **Solution** the NV-Measurement



The prepared fluid enters the NV-Measurement chip from the port on the top left and is then routed over the top of the NV-Diamond. The Sample is analysed by the NV-Setup and then discarded.

This glass chip also gets cleaned automatically leaving no need for human labour.

# Our Solution the Chip -Stack




This view shows how the two chips fit together inside the actual machine. The laser for pumping the diamond and the emitted light is the basis of the analysis are also visible.







# Disruptive new Diagnoses

- Diagnoses of Alzheimers
- Effective monitoring of chronic conditions
- Early detection of cancer
- Beter dosage of medicine







# Our Competitors (diagnoses of Alzheimers)

	 BioQube	 SIEMENS Healthineers	 Roche	 GE HealthCare
Cheap	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Fast	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Non-Invasive	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Accurate	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Early diagnoses	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>
No need for expertly trained staff	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

# Market Size




# Roadmap



Helping 33 million people with Alzheimers Disease

Learn more:



Tobias Beck  
Hardware Design



Matthäus Hirsch  
Hardware Design



Tim Maehrholz  
Market Research



Lisa Haag  
Biological Research



Jascha Fricker  
Technical Design



Benedict Brouwer  
Technical Design